

MURLIE-X

Installation Guide

© 2026 NVISAGE All rights reserved.

Under copyright laws, this manual should not be reproduced, in whole or in part, without the written permission of NVISAGE. The specification or any information can change without notice.

NVISAGE
TECH
NOLO
GIES

Table of Contents

- 1. Instruction 3
 - 1.1 Instruction 4
 - 1.1.1 Target Audience 4
 - 1.1.1 Symbols 4
 - 1.1.2 Revision History 4
 - 1.1.3 Reference Documents 5
- 2. Installation 6
 - 2.1 Installation Method 7
 - 2.2 FAQ 9

1 Instruction

NVISA
GE
TECH
NOLO
GIES

1.1 Instruction

1.1.1 Target Audience

This document is intended for users who need the information on MURIE-X installation.

1.1.1 Symbols

This product should be operated under the safety instructions with the warning or caution symbol in this manual. It is important for you to read and understand the contents to operate the products safely.

Caution



- This symbol is used for statements intended to prevent potentially hazardous situations that may result in death, personal injury, or significant property damage. Users should fully understand the information associated with this symbol..

Information



- This symbol is used to indicate product-related notes and supplementary information. Users are advised to read statements containing this symbol carefully.

1.1.2 Revision History

Version	Date	Description
1.0	2026-06-17	Initial release

1.1.3 Reference Documents


The QR codes and URLs related to MURLIE-X are as follows:

Document	QR Code	URL
NEO-E Quick Start Guide		https://qrfy.io/p/qBX0iWScNI
NEO-E User Manual		https://qrfy.io/r/zzk7HHtogP
N-Bridge User Manual		https://qrfy.io/r/LK9T_uBNNA
N-Bridge OneStop Settings Guide		https://qrfy.io/r/i0QqmNu_BM
N1 User Manual		https://qrfy.io/r/cR4VSiQIKy

NVISAGE TECH NOLO GIES



2.1 Installation Method

- 1 When you run the **Install.bat** program, a console window will appear and the software will be installed automatically. Once the installation is complete, the message "Installation Complete!" will be displayed in the console window.




- 2 The program is installed in the directory **C:\Program Files\NVisage**, and open the **Murlie-X** file on the desktop. If an error occurs indicating that the DirectX program is not installed, download it from the following link:

<https://www.microsoft.com/ko-kr/download/details.aspx?id=35>



- 3 After running MURLIE-X, connect the sensor with N-Bridge version 1.5.1 or higher, then click the Authentication button on the screen below.



- 4 Fill in the required fields below and click the **Request Code** button.



- First Name and Last Name can each contain up to 20 characters.
- If the same user performs authentication more than twice at a commercial site, it is recommended to use a brand name or company name instead of a personal name.

5 An authentication email will be sent to the email address. Enter the received code in the Authentication Code field below and click Submit. If you do not receive the email, please check your spam folder.

6 Once the code is verified, the message "The authentication code has been successfully verified" will appear, and the Register Product button will be displayed. Click the button to complete the registration process.

2.2 FAQ

Q: Can multiple sensors be registered to one PC?

A: No. Only one sensor can be registered per PC. If you change the sensor, you must re-authenticate and register the new sensor.

Q: Can one sensor be registered to multiple PCs?

A: Yes. If the PCs are different, a single sensor can be registered to multiple PCs.

Q: What happens if a different sensor is registered on a PC where a sensor is already registered?

A: When a new sensor is registered, the previously registered sensor information is automatically deleted and replaced with the new sensor information.

Q: What should I do if I replace the sensor?

A: If the sensor is changed, the existing registration cannot be used. You must re-authenticate and register the new sensor.

